

As always.
Only better.

The counter-running Pütter bandaging technique causes very high decongestion which also reaches deep veins.

Tips for use

- After washing (up to 95°C) the bandage regains its full elasticity.
- More information including an application video available at www.puetterflex.com.

How to apply a "modified Pütter bandage"



- 1** Place the foot at a right angle to the leg and start the first turn of the bandage at the metatarsophalangeal joints and run outwards.
- 2** After 2-3 turns around the mid-foot, the bandage covers the heel and comes back to the instep over the medial ankle.
- 3** With each turn, first above and then below, the edges of the first heel wrap are additionally secured. After a further turn around the midfoot, the bandage is taken back to the ankle and across the ankle joint line.
- 4** The bandage then conforms to the leg spirally, widely spread, in order to enclose the calf, drawn only in the direction of application.



- 5** From the popliteal space, the bandage is wrapped once around the leg and from there follows the shape of the leg back down, covering any gaps in the bandage.
- 6** The second bandage is applied in the opposite direction starting from the outside towards the inside of the ankle around the heel.
- 7** The further two turns fix first the upper and then the lower edge of the heel wrap.
- 8** Finally, the bandage is wrapped once more around the midfoot and then spiralled upwards and back down again in the same way as the first one. The completed bandage is fixed with adhesive tape (e.g. Omniplast) on the outside of the bandage.



PAUL HARTMANN AG
D-89522 Heidenheim
Germany

Visit our Website
www.hartmann.info



Bandages in various widths for individual use as needed.

Presentations	Unit	Code No.	Units per Case
PütterFlex with bandage clips, stretched length 5 m, individually boxed			
6 cm	1 bandage	931827	40
8 cm	1 bandage	931828	40
10 cm	1 bandage	931829	40
12 cm	1 bandage	931830	60
PütterFlex with bandage clips, stretched length 5 m, loose in boxes			
6 cm	10 bandages	931831	8
8 cm	10 bandages	931832	8
10 cm	10 bandages	931833	8
12 cm	10 bandages	931834	8
PütterFlex-set			
10 cm	2 bandages	931826	30



NEW!

Pütter®Flex
Comfort compression for more mobility.



Compression Therapy

Compression and comfort – No longer a contradiction.

In order to be successful, compression treatment requires the active cooperation of the patient. Therefore the ideal bandage is one which is comfortable to wear and ensures efficient compression at the same time. These properties substantially increase the chances that the patient will continue treatment until the end.

This combination of comfort and compression was, however, not always available in the past. Especially when the foot was wrapped with bandages, it often did not fit into the shoe anymore and inhibited the patient's mobility. In addition, the pressure necessary for treatment was often found to be uncomfortable. Now there is a solution to these problems: the new PütterFlex.



Healthy mobility

Since the ankle joint remains as mobile as possible, the patient is still able to move about freely in everyday life. In addition, the efficacy of compression treatment is enhanced because the muscular venous pump is supported by every step.

The air-permeable material promotes a comfortable feeling in wear for the patient.



Pütter®Flex

The bandage can also be worn at night due to its low resting pressure.

Since the bandage is flexible in all directions, it adapts better even to difficult areas. This makes application easier and increases wearing comfort.



A perfect fit!

Because the bandage is very thin, the wrapped foot fits better in a shoe. As a result, the treatment is better suited for daily use.

Benefits for the patient

- High wearing comfort
- Fits in normal shoes
- More mobility in the ankle joint
- Can be worn day and night

Effectively combined. Better results.

For example PütterFlex with a HARTMANN wound dressing individually adapted to the healing phase. This highly effective treatment can be specifically adapted from case to case and makes successful treatment easier.



Products. Knowledge. Experience.
Combine for maximum benefit.

-  **TenderWet 24 active**
For active cleansing of venous leg ulcers.
-  **Atrauman Ag**
For treatment of wounds in danger of infection or infected wounds as well as for prophylaxis against infections in chronic or acute wounds.
-  **PermaFoam**
For treatment of heavily to moderately exuding leg ulcers in the cleansing phase and infected wounds in combination with Atrauman Ag.
-  **Hydrocoll**
For treatment of moderately to lightly exuding, non-infected wounds.
-  **Hydrocoll thin**
For treatment of lightly exuding and already epithelizing wounds.
-  **Hydrosorb**
For treatment of ulcers in the epithelization stage.
-  **HydroTac**
For treatment of moderately to lightly exuding ulcers.

www.combisensation.com